

Experience high power polymerizations  
and whitening treatments



# *Ti-Lite*

**Cordless LED Curing Light**

# GT-1500

**MONITEX**

# Ti-Lite GT-1500

## Cordless LED Curing Light

Ti-Lite GT-1500 is a powerful curing light with a light intensity of 1,500 mW/cm<sup>2</sup>, ensuring quick and complete polymerization of dental composites and resins. In addition, GT-1500 can also be used in conjunction with whitening gels to perform whitening treatments, making it a handy unit to have in the clinic.

### 3 Curing Modes

#### 1 | Full Power Mode

Full power for optimal high speed polymerization.

**Cycle Duration** Default 8 secs (adjustable)



#### 2 | Whitening Mode

15 minute whitening cure at 70% power.


**Cycle Duration** Default 15 mins



#### 3 | Pulse Mode

Unit turns on and off every second to reduce heat production.

**Cycle Duration** Default 20 secs (adjustable)



**Bleaching Arch Included**

#### Specifications:

- LED Lamp
  - Light Intensity
  - Wavelength Range
  - Fiber-optic Probe
  - Handpiece Dimensions
  - Handpiece Weight
  - Charger Stand Dimensions
  - Charger Stand Weight
  - Power Input
  - Power Output
- 5 W  
1,500 mW/cm<sup>2</sup> and up  
420~490 nm (peaks at 455~465 nm)  
Ø10~Ø8mm, autoclavable up to 134°C  
185(L) x 38(W) x 82(H) mm  
206 g  
193(L) x 83(W) x 95(H) mm  
177 g  
AC 100~240V, 50/60 Hz  
DC 9V  $\approx$  1.3A



**MONITEX** INDUSTRIAL CO., LTD

5F, 70, Sec. 1, Guangfu Rd., Sanchong Dist., New Taipei City 24158, TAIWAN  
TEL: + 886 2 29993869 FAX: + 886 2 29993647  
Email: [monitex@monitex.com.tw](mailto:monitex@monitex.com.tw)  
[www.monitex.com.tw](http://www.monitex.com.tw)

